



EZ-Vision™

A Fluorescent Dye for the Instant Visualization of DNA Bands in Agarose Gels

Supplied in 6X Loading Buffer

| <u>Code</u> | <u>Description</u> | <u>Size</u> |
|---------------|---|-------------------------------------|
| N313-KIT | EZ-Vision™ Three, DNA Dye as Loading Buffer, 6X Supplied in Loading Buffer, 6X containing 15% Ficoll Includes 3 tracking dyes migrating at 10, 400 and 4000 bp | 5 x 1.0 ml |
| N313-Q-SAMPLE | EZ-Vision™ Three, DNA Dye as Loading Buffer, 6X Supplied in Loading Buffer, 6X containing 15% Ficoll Includes 3 tracking dyes migrating at 10, 400 and 4000 bp | 0.3 ml |
| N472-KIT | EZ-Vision™ One, DNA Dye as Loading Buffer, 6X Supplied in Loading Buffer, 6X containing 15% Ficoll Includes a single tracking dye migrating at 10 bp | 5 x 1.0 ml |
| N472-Q-0.5ML | EZ-Vision™ One, DNA Dye as Loading Buffer, 6X Supplied in Loading Buffer, 6X containing 15% Ficoll Includes a single tracking dye migrating at 10 bp | 0.5 ml |
| N473-2PK | EZ-Vision™ Sample Kit Includes: EZ-Vision™ Three, DNA Dye as Loading Buffer, 6X EZ-Vision™ One, DNA Dye as Loading Buffer, 6X | 1 Kit 1.0 ml 1.0 ml |
| N473-3PK | EZ-Vision™ Sample Kit Includes: EZ-Vision™ One, DNA Dye as Loading Buffer, 6X EZ-Vision™ two, DNA Dye as Loading Buffer, 6X EZ-Vision™ Three, DNA Dye as Loading Buffer, 6X | 1 Kit 1.0 ml 1.0 ml 1.0 ml |
| N650-KIT | EZ-Vision™ Two, DNA Dye as Loading Buffer, 6X Includes 2 tracking dyes migrating at 400 and 4000 bp | 5 x 1.0 ml |
| N650-Q-SAMPLE | EZ-Vision™ Two, DNA Dye as Loading Buffer, 6X Includes 2 tracking dyes migrating at 400 and 4000 bp | 0.3 ml |

General Information:

EZ-Vision™ is a non-mutagenic fluorescent reagent that produces instant visualization of DNA bands upon UV illumination of agarose gels. EZ-Vision™ is non-mutagenic with no hazardous shipping, handling or disposal costs. Supplied in AMRESCO's 6X DNA Loading Buffer, EZ-Vision™ forms a tight complex with the sample DNA and co-migrates with it during electrophoresis. Post-run staining and destaining is completely eliminated and results can be visualized immediately after the run by placing the gel on a standard UV transilluminator. It is ideal for applications needing rapid DNA band visualization and for environments requiring a safe, non-hazardous alternative to Ethidium Bromide.

EZ-Vision™ is available in 3 versions that differ only by the tracking dyes included in the loading buffer. EZ-Vision™ Three contains 3 tracking dyes that migrate at 4,000 bp, 400 bp, and 10 bp. EZ-Vision™ Two contains 2 tracking dyes that migrate at 4,000 bp and 400 bp. EZ-Vision™ One contains only a single fast-running tracking dye that migrates at approximately 10 bp in a 1% agarose gel.



Storage/Stability:

EZ-Vision™ is stable for at least 1 year at 2 - 8°C.

EZ-Vision™ is light sensitive and should be stored protected from light. Normal usage can be carried out under ambient light.

Spectral Information:

Excitation = 364 nm

Emission = 454 nm

Application Disclaimer

For Research Use Only. Not for Therapeutic or Diagnostic Use.

EZ-Vision™ tracking dyes:

| EZ-Vision™ Type | Color | Mobility (1% Agarose) |
|------------------------|--------------|----------------------------------|
| EZ-Vision™Three | Light blue | 4000 bp |
| | Dark blue | ~400 bp |
| | Magenta | ~ 10 bp |
| EZ-Vision™Two | Light blue | 4000 bp |
| | Dark blue | ~400 bp |
| EZ-Vision™One | Magenta | ~ 10 bp |



Protocol:

No mutagenic or genotoxic effects are observed in AMES testing or sister-chromatid exchange assays of EZ-Vision™. In addition, EZ-Vision™ passes the CCR Title 22 Fathead Minnow Hazardous Waste Screen Bioassay. EZ-Vision™ is not toxic at the concentrations used but standard handling precautions are advised for all nucleic acid binding reagents. All local regulations should be followed when using and disposing of this reagent.

1. Vortex EZ-Vision™ for 30 seconds prior to use.
2. Dilute 1 part EZ-Vision™ with 5 parts DNA sample and mix.
 - Note: EZ-Vision™ must be added to DNA markers in order to visualize the ladder bands simultaneously with the sample after electrophoresis.
3. Load sample and run according to standard procedure.
4. After the run, remove gel and place on UV transilluminator to immediately visualize bands. DNA bands will emit a whitish-blue fluorescence against a dark background. Fluorescence should be visible on a transilluminator for at least 24 hours after electrophoresis if the gel fluorescence has not been bleached.
5. Gels can be post stained with Ethidium Bromide if desired.

Gel Documentation

- **Black and White Polaroid Photography:** EZ-Vision™ stained gels can be photographed with a standard UV transilluminator, Polaroid #667 film and filters used to photograph green dyes such as SyBR™ Green. Longer exposure times are required for EZ-Vision™ stained gels since Polaroid #667 film is not optimized for sensitivity in the blue emission range. Filters typically used to photograph Ethidium Bromide (EtBr) stained gels can be used, but exposure times should be doubled or tripled to obtain sufficient contrast to represent the image that is visually perceived.
- **Gel Imaging Systems:** EZ-Vision™ stained gels are compatible with digital imaging systems. Please contact your system manufacturer with the excitation and emission information listed below to obtain information on appropriate filters.

Excitation = 364 nm

Emission = 454 nm

Downstream Applications

DNA stained with EZ-Vision™ is compatible with a variety of downstream applications including ligation reactions, transformation procedures and PCR amplification.

Results Summary

Ligation, transformation and PCR amplification procedures were tested on parallel samples of plasmid DNA fragments stained with either EZ-Vision™ or EtBr.

- **Recovery from gel slices:** The EZ-Vision™ reagent does not interfere with recovery of DNA fragments from agarose gels. Comparison studies from 1% agarose gels were performed to determine the recovery of EZ-Vision™ stained DNA fragments versus EtBr stained fragments. In these tests, a DNA aliquot was combined with EZ-Vision™ and applied to a 1% TAE agarose gel. A second aliquot was loaded on a 1% TAE agarose gel containing 50 µg/ml EtBr. Gels were run according to standard procedures. An 850 bp DNA band was excised from each gel and purified on QIAquick® Gel Extraction Kits (Qiagen Group, Hilden, Germany). Equal aliquots were applied to 1% TAE agarose gels. DNA recoveries were determined by band intensity after electrophoresis. There was no significant difference in the amount of DNA recovered from each gel.
- **Ligation and Transformation Efficiency:** Ligation and transformation efficiency of the EZ-Vision™ stained 850 bp fragment is similar to the EtBr stained fragment. In parallel reactions, the purified 850 bp fragments were ligated into an Ampicillin resistant vector and transformed into E.coli. The transformed cultures were plated onto LB-Ampicillin media. The total number of colonies for each sample was determined after overnight incubation. The number of Ampicillin resistant colonies was similar for each sample.
- **PCR Amplification:** Amplification of an EZ-Vision™ stained DNA fragment was equivalent to the EtBr stained fragment. A 1.3 kb DNA fragment stained with either EZ-Vision™ or EtBr was amplified by PCR using Pfu polymerase. The products were applied to a 1% TAE agarose gel. The 1.3 kb band was excised and purified as described above (see Recovery from gel slices). The recovered DNA samples were quantitated by a NanoDrop® Spectrophotometer (NanoDrop® Technologies, Wilmington, DE). The amount of DNA recovered from each sample was equivalent indicating that EZ-Vision™ stained DNA was a suitable template for PCR reactions.
- **Sequencing:** The use of EZ-Vision™ stained DNA as a sequencing template is currently being evaluated.



Related Products

| Code | Product |
|----------------|---|
| Agarose | |
| 0710-500G | Agarose I™, 500 g General Use (also available as tablets, K857-100TABS) |
| J234-250G | Agarose SFR™ Super Fine Resolution for superior resolution of nucleic acid fragments between 200 and 1,000 base pairs. |
| E776-100G | Agarose 3:1 HRB™ High Resolution Blend for resolution of nucleic acid fragments below 1,000 base pairs. |
| Buffers | |
| 0658-4L | TBE Buffer, 10X Liquid Concentrate |
| 0478-2PK | TBE Buffer, 10X Ready-Pack™ |
| 0796-1.6L | TAE Buffer, 25X Liquid Concentrate |
| Markers | |
| K180-250UL | 100 bp Ladder 13 bands ranging from 100-3000 base pairs |
| E854-50RXN | PCR DNA Marker™ 8 bands ranging from 50 to 2000 base pairs |
| K181-500UL | 1 kb Ladder 11 bands ranging from 500 to 10,000 base pairs |

Visit the AMRESCO website to view
additional related products

www.amresco-inc.com



www.interchim.com



211 bis Avenue Kennedy - BP 1140
03103 Montluçon - France
33 (0) 4 70 03 88 55
Fax 33 (0) 4 70 03 82 60
e-mail interchim@interchim.com

Agence Paris - Normandie
33 (0) 1 41 32 34 40
Fax 33 (0) 1 47 91 23 90
e-mail interchim.paris@interchim.com

References:

1. Andrews, A.T. Electrophoresis: Theory, Techniques, and Biochemical and Clinical Applications 2nd ed., New York, (1988), 21-24.
2. Ogden, R.C. and Adams, D.A. Electrophoresis in agarose and acrylamide gel. Methods Enzymol., 152, 61-87 (1987)

ZY0390 Rev.7 9/20/2009 ©Copyright 2006, 2007, 2008, 2009 by AMRESCO Inc.

All Rights Reserved.

AMRESCO® is a registered trademark of AMRESCO Inc.

Qiaquick® is a registered trademark of the QIAGEN Group.

Tiffen® is a registered trademark of The Tiffen Company.

Nanodrop® is a registered trademark of Nanodrop® Technologies.

