

Product Specification Sheet

Product Name	oYo-Link [®] mIgG1 VcMMAE
Catalog #	AT7001-mIgG1
Lot #	Lot # provided when product shipped
Molecular Weight	10 kDa
Concentration	337 µg/mL
Molar Concentration	33 µM
Storage Buffer	1 x PBS
Shipping Condition	This product is shipped as a dried clear or white pellet at ambient temperature.
Reconstitution	Before opening, briefly centrifuge each tube to ensure that the entire pellet is at the bottom. Add 100 µL sterile ddH₂O and vortex tube briefly until fully dissolved.
Storage and Stability	The dried product can be stored at 4°C or -20°C upon receipt. Following reconstitution, the product can be stored at 4°C for short-term durations (< 1 week) and -20°C to -80°C for longer time-spans.
Conjugation Instructions	After reconstitution, each tube is capable of labeling 100 µg of only Mouse IgG1 antibody . Please find and follow the Antibody Conjugation User Manual at www.alphathera.com/user-manuals for conjugation instructions.
Notes	<p>Disposal of product should be performed in compliance with all applicable federal, state, and local regulations.</p> <p>This product does not meet hazard classification criteria based on evaluations made by our company under the OSHA Hazard Communication Standard, 29 CFR 1910.1200.</p>

For more information, please visit www.alphathera.com or contact support@alphathera.com.